

Via A. Volta n. 27 Tel. 39 02 90093082 www.gambetti.it

20082 Binasco (Milano) Fax. 39 02 9052778 info@gambetti.it www.plasmi.eu





PORTABLE VACUUM CALIBRATION SYSTEM

The PVS6E consists of up to three MKS Baratron® Type 690 High Accuracy Absolute capacitance manometers as transfer standards and a high vacuum pumping system, mounted in a sturdy transportable cart. The PVS6E provides NIST traceable calibrations over the pressure range from 10⁻⁵ to 1000 mmHg and can be used to calibrate MKS Baratron® Capacitance Manometers, thermocouple gauges, Pirani gauges, convection enhanced Pirani gauges, other capacitance manometers and transmitters, and mechanical/ dial gauges. Hot or cold cathode ionization gauges can be calibrated over the upper end of their range.

The PVS6E is a completely re-engineered version of the popular PVS6D. Standard features include a narrow chassis for ease of movement through confined areas, clean room compatibility, dry pumping system, and two-channel display/readout.

Features & Benefits

- Wide 10⁻⁵ to 1000 mmHg calibration range for many different types of vacuum gauges
- · Calibrates instruments in situ to keep process up and running and to reduce maintenance
- Simple, easy to use system encourages regular instrument calibration resulting in higher process quality and yields
- Uninterruptible power supply keeps standards ready to use

- · Can be used as a process or instrument diagnostic tool to eliminate unnecessary replacement of functioning instruments
- Manual operation to meet your requirements
- · NIST traceable calibration helps comply with government regulations, ISO 9000 and other QA program requirements





The PVS6E pumping system consists of an air-cooled hybrid turbomolecular pump, backed by a two stage diaphragm pump. A high vacuum gauge is provided to monitor system performance and base pressure to ensure pressure is low enough for proper zeroing of the transfer standards. The PVS6E is available with manual pressure control capability. The PVS6E has an uninterruptible power supply with sufficient storage to keep the transfer standards warmed up and ready to use for up to 60 minutes, allowing ample time to transport the PVS6E to the process system.

The PVS6E can perform off-line as well as *in situ* calibrations (see Figures 1 and 2).




Figure 1 — Configuration for off-line calibration




Figure 2 — Configuration for in situ calibration with the process gauge mounted on the tool's chamber



Specifications

Transfer Standard MKS Type 690 Absolute Sensor (NIST traceable); see 690 data sheet for other specifications

Baratron® Pressure Sensor Accuracy

(non-linearity, hysteresis, and non-repeatability)

0.08% of Reading (± temp. coefficients), 0.05% of Reading (± temp. coefficients)

or optional for 1-1000 range

Resolution (of F.S.) 1 x 10⁻⁶

Unit Port Fittings8 VCR® (3 each)Operating Temperature Range15° to 40°C

Readout Electronics MKS Type 670, 5½-digit LCD

Power Required 110 VAC, 60 Hz; 220 VAC, 50 Hz; 15 Amps

Warm-up Time/Start Time 4 hours. However, use of the internal 500 VA UPS will keep the Transfer Standard warmed up

and ready for immediate use for at least 60 minutes.

Pressure Control Manual needle valve to control gas inlet; manual isolation valve and metering bypass

valve to control pump speed

Vacuum Pumping System Mechanical pump: 18 Liters/min diaphragm

High vacuum pump: hybrid turbo/drag air-cooled

Normal base pressure: <5 x 10⁻⁷ mmHg

Vacuum Pump Discharge Fitting NW 16 KF

High Vacuum Gauge IMAG ionization gauge and controller provides 1 x 10⁻³ to 1 x 10⁻⁹ mmHg

measurement

 Test/Calibration Gas
 5-10 psig dry filtered N₂

 Test Gas Inlet Fitting
 ¼" Swagelok® bulkhead

 Test Port Fitting (3)
 Swagelok® 8 VCR® male (3)


Standby Power 500 V.A. uninterruptible power supply keeps three Transfer Standards warmed up and ready to


use for up to 60 min.

Weight Approximately 250 lbs (110 kg)



Ordering Information







Dimensional Drawing —

Note: Unless otherwise specified, dimensions are nominal values in inches (mm referenced).



Ordering Code Example: PVS6EXXXYZ

Type PVS6E Portable Vacuum Calibration System	Code	Configuration
Type Number	PVS6E	PVS6E
Flow Scale Pressure Ranges in mmHG (XXX)		
0.1 1 10 100 1000 Vacant channel	A B C D E	CDO
Input Power (Y)		
Standard USA & Canada input power 110 VAC, 60 Hz 220 VAC, 50 Hz	U F	U
Accuracy (Z)		
Best Available ±0.08% of Reading	A B	А

MKS Instruments, Inc. Global Headquarters

2 Tech Drive, Suite 201 Andover, MA 01810

Tel: 978.645.5500 Tel: 800.227.8766 (in U.S.A.)

Tel: 800.227.8766 (in U.S Web: www.mksinst.com

MKS Instruments, Inc. Pressure & Vacuum Measurement Solutions

Six Shattuck Road Andover, MA 01810 Tel: 978.975.2350

Some Baratron® capacitance manometer products may not be exported to many end user countries without both US and local government export licenses under ECCN 2B230. MKS products provided subject to the US and other country export regulations. Diversion or transfer contrary to U.S. or other export laws is prohibited.

Specifications are subject to change without notice. mksinst™ is a trademark and Baratron® is a registered trademark of MKS Instruments, Inc., Andover, MA. Swagelok® and VCR® are registered trademarks of Swagelok Marketing Co., Solon, OH. DeviceNet™ is a trademark of the Open DeviceNet Vendor Association, Coral Springs, FL.